

Technical Data

Product: DIVAC 1.2 L
Part No.: 13508

Inlet connection:	Hose nozzle ID 10, G 1/4" Inside thread	
Outlet connection:	Hose nozzle ID 10, G 1/4" Inside thread	
Pumping speed:	1.2 m ³ /h	0.7 cfm
Ultimate pressure:	≤ 8 mbar	≤ 6 Torr
Max. outlet back pressure: (absolute)	2000 mbar	1500 Torr
Noise level:	50 dB(A)	
Mains connection:	1-ph, 100 V, ±10%, 50/60 Hz	
Motor rating:	120 W	0.17 HP
Rotation speed at 50 Hz:	1500 min ⁻¹	1500 rpm
Rotation speed at 60 Hz:	1800 min ⁻¹	1800 rpm
Motor protection rating:	IP 44	
Admissible gas inlet temperature:	+5 to +40°C	+41 to +104°F
Admissible ambient temperature:	+5 to +40°C	+41 to +104°F
Dimensions:	see dimension sheet	
Weight:	9.3 kg	20.5 lbs

Technical data are subject to change